
Name	Class




1 Label the microscope using words from the box.



eyepiece lens
focusing wheel
mirror
objective lens
slide
stage

2 Use these clues to complete the crossword.



Across

- **4** Lens closest to the slide.
- **5** The specimen is placed on this.
- **6** You need this to see something

Down

- 1 Lens that you look down.
- **2** What the slide is placed on.
- **3** A wheel to make the image clear.
- **3** Here are some instructions on using a microscope. Put the number '1' next to the instruction you should do first. Put the number '2' next to the instruction you should do next and so on.

	Look	into	the	eyepiece	lens.
--	------	------	-----	----------	-------

_____ Place the smallest objective lens over the hole in the stage.

_____ Place the slide on the stage.

_ Turn the coarse focusing wheel until what you see is clear.

Turn the coarse focusing wheel to make the objective lens as

close to the stage as possible.

____ Adjust the light source.

